



Medline International

www.medline.com/uk

Manufactured by



Technical Sheet

Triflex[®] LP

Sterile Latex Powdered Surgical Gloves by Cardinal Health

Reference(s), size(s)

2D7255LP	2D7260LP	2D7265LP	2D7270LP	2D7275LP	2D7280LP	2D7285LP	2D7290LP
5,5	6	6,5	7	7,5	8	8,5	9

Primary material

Natural rubber latex

Powder content: <15.0 milligrams per square decimeter, in accordance with ASTM D 6124 and ISO 21171.

Caution: This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses.

Color

Cream

Chemical additives

ZDBC (Zinc Dibutyldithiocarbamate) is used in manufacturing and below detectable level.

Donning agent

Absorbable dusting powder

Leachable protein

<200 micrograms/gram, in accordance with EN 455-3 using ASTM D 5712 (Modified Lowry Method).

Caution: Safe use of these gloves by or on latex-sensitive individuals has not been established.

Thickness
(per ASTM D 3577 \geq 0,10 mm)

Finger Tip	0,22 mm ⁽¹⁾
Palm	0,20 mm ⁽¹⁾
Cuff	0,16 mm ⁽¹⁾

Length (per EN 455 > 270 mm)

305 mm^(x)

Force at break, before accelerated ageing (per EN 455-2 \geq 9 N)

18 N

Force at break, after accelerated ageing (per EN 455-2 \geq 9 N)
(7 days at a temperature of 70°C in an oven)

19 N

Elongation at break, before accelerated ageing
(per ASTM D 3577 \geq 750 %)

873 %

Elongation at break, after accelerated ageing (per ASTM D 3577 \geq 560 %)
(7 days at a temperature of 70°C in an oven)

883 %

Freedom from holes
(per EN 455 AQL 1.5)

0.65 AQL

Viral penetration

Tested and passed, in accordance with ASTM F 1671



STERILE R



Medline International

www.medline.com/uk

Manufactured by



Technical Sheet

Triflex[®] LP

Sterile Latex Powdered Surgical Gloves by Cardinal Health

Sterilization	Gamma radiation
Design	Mould has an independent thumb for more natural movement in the fingers, thumb and palm.
Cuff design	Plain (Non Beaded)
Expiration date	5 years from date of manufacture. Manufacture and expiration dates are printed on packaging (YYYY-MM format).
Packaging	50 pairs per dispenser box / 4 boxes per carton
Storage recommendations	Protect from freezing. Avoid excessive heat. Keep dry. Product should be shielded from direct sunlight, fluorescent lighting, X-rays, moisture and ozone. Do not store in temperatures above 40°C.
Country of origin	Thailand
Regulations and Quality standards	Cardinal Health manufacturing facilities are certified to EN ISO 13485 by BSI. Products conform to EU Directive 93/42/EEC and CE marked. Products meet the requirements of European harmonized standards EN 455-1, -2, -3 and -4.

⁽¹⁾ Test results by Cardinal Health (per ASTM D 3577)

^(X) Data average was calculated using size 7½

