



-  5 mm Trode
-  20 mm Trode
-  35 mm Trode
-  80 mm Trode

**GENERAL RECOMMENDATIONS:**  
 1-3 treatments;  
 delivered 2-3 weeks apart

• **Back**  
 • **Tendon/Ligaments**  
 1-3 treatments;  
 delivered 4 weeks apart

**OTHER:**  
 • **Splint**  
 5mm/20mm  
 1500 pulses  
 E6

• **Superficial Wounds**  
 5mm  
 8 pulses/cm<sup>2</sup>; min 500  
 E1-E2